

## Medical Assembly & Test Machine Examples:

**Catheter Assembly & Test**– Feed & cut tubing to various lengths, die cut eyes, tip tube end, apply CA adhesive and add fitting. Test flow, vision inspect, offload. 100% part inspection, 4– styles of parts with various tube lengths and diameters. 4– second cycle. Intermittent motion dial with cam lift tooling.



**Oral Hygiene Appliance Assembly & Test (Suction Toothbrush)**– Feed & assemble brush end wand to fitting with CA adhesive, install foam pad with adhesive, check flow and assembly. 100% part inspection, 2– styles of parts, 2.5 second cycle per part on a precision-link conveyor.

